Estudio **ANGIOLITE**

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Comentarios del experto

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First-in-man randomised comparison of the Angiolite durable fluoroacrylate polymer-based sirolimus-eluting stent versus a durable fluoropolymer-based everolimus-eluting stent in patients with coronary artery disease: the ANGIOLITE trial

EuroIntervention. 2019 Dec 20;15(12).

AIMS

The durable fluoroacrylate polymer-based sirolimus-eluting stent (Angiolite SES) has shown promising preclinical and clinical results regarding inflammatory vascular reaction and neointimal healing. We aimed to compare performance between the Angiolite SES and an everolimus-eluting stent (EES) in patients with coronary artery disease.

METHODS AND RESULTS

The ANGIOLITE trial, a prospective, randomised, multicentre trial, compared the restenosis parameters of both stents in de novo coronary lesions. The primary endpoint was late lumen loss at six-month angiographic follow-up. In-stent healing was assessed by optical coherence tomography (OCT). The main clinical endpoint was target lesion failure (TLF) evaluated up to 24 months.

A total of 223 patients were randomised 1:1 to EES or SES. At six months, in-stent late lumen loss was 0.08 mm (\pm 0.38) for EES vs 0.04 mm (\pm 0.39) for SES (difference = -0.04 mm, 95% CI: -0.15, 0.07, p for non-inferiority = 0.002). By OCT, the rate of uncovered to total number of struts score > 30% was comparable between the groups whereas neointimal thickness was reduced in the SES arm (9.0% [7.6, 10.6] vs 9.9% [8.5, 11.3], p = 0.41; and 86.4 [81.6, 91.2] μ m vs 72.1 [68.2, 76.0] μ m, p < 0.01, respectively).

At 24 months, TLF occurred in eight patients (7.6% [3.3, 14.5]) in the EES arm and in seven patients (7.1% [2.9, 14.0]) in the SES arm (p = 0.88). The definite/probable stent thrombosis rate was comparable between the groups (1.9% [0.2, 6.7] vs 1.0% [0.0, 5.5] EES vs SES, respectively; p = 0.59).

CONCLUSIONS

This trial demonstrates similar antirestenotic efficacy at midterm follow-up of the Angiolite SES vs an EES. Clinical endpoints were comparable between the groups at two-year follow-up.

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
1. Hartley A, Marshall DC, Saliccioli JD, Sikkil MB, Maruthappu M, Shalhoub J. Trends in mortality from ischemic heart disease and cerebrovascular disease in Europe: 1980 to 2009. *Circulation*. 2016;133(20):1916-26.
2. Ibáñez B, James S, Agewall S, Antunes MJ, Bucciarelli-Ducci Ch, Bueno H, et al. 2017 ESC Guidelines for the management of acute myocardial infarction in patients presenting with ST-segment elevation: The Task Force for the management of acute myocardial infarction in patients presenting with ST-segment elevation of the European Society of Cardiology (ESC). *Eur Heart J*. 2018;39(2):119-77.
3. Brophy JM, Belisle P, Joseph L. Evidence for use of coronary stents. A hierarchical bayesian meta-analysis. *Ann Intern Med*. 2003;138:777-86.
4. Joner M, Finn AV, Farb A, Mont EK, Kolodgie FD, Ladich E, et al. Pathology of drug-eluting stents in humans. Delayed healing and late thrombotic risk. *J Am Coll Cardiol*. 2006;48(1):193-202.
5. Tada T, Byrne RA, Simunovic I, King LA, Cassese S, Joner M, et al. Risk of stent thrombosis among bare-metal stents, first-generation drug-eluting stents, and second-generation drug-eluting stents: Results from a registry of 18,334 patients. *JACC Cardiovasc Interv*. 2013;6:1267-74.
6. Smits PC, Hofma S, Togni M, Vázquez N, Valdes M, Voudris V, et al. Abluminal biodegradable polymer biolimus-eluting stent versus durable polymer everolimus-eluting stent (COMPARE II): A randomised, controlled, noninferiority trial. *Lancet*. 2013;381:651-60.
7. Bona KH, Mannsverk J, Wiseth R, Aaberge L, Myreng Y, Nygard O, et al; NORSTENT Investigators. Drug-eluting or bare-metal stents for coronary artery disease. *N Engl J Med*. 2016;375:1242-52.
8. Moreu J, Moreno-Gómez R, Pérez de Prado A, García del Blanco B, Trillo R, Pinar E, et al. First-in-man randomized comparison of the Angiolite® durable fluoroacrylate polymer-based sirolimus-eluting stent versus durable fluoropolymer based everolimus eluting stent in patients with coronary artery disease: The ANGIOLITE trial. *EuroIntervention*. 2019;15(12):e1081-e1089.
9. Joner M, Nakazawa G, Finn AV, Quee SC, Coleman L, Acampado E, et al. Endothelial cell recovery between comparator polymer-based drug-eluting stents. *J Am Coll Cardiol*. 2008;52(5):333-42.
10. Spaulding C, Teiger E, Commeau P, Varenne O, Bramucci E, Slama M, et al. Four-year follow-up of TYPHOON (trial to assess the use of the CYPHer sirolimus-eluting coronary stent in acute myocardial infarction treated with Balloon angioplasty). *JACC Cardiovasc Interv*. 2011;4(1):14-23.
11. Estévez-Loureiro R, Pérez de Prado A, Pérez-Martínez C, Cuellas-Ramón C, Regueiro-Purriños M, Gonzalo-Orden JM, et al. Safety and efficacy of new sirolimus-eluting stent models in a preclinical study. *Rev Esp Cardiol (Engl Ed)*. 2015;68(12):1118-24.
12. Puri R, Otaegui I, Sabaté M, Serra-Peñaranda A, Puigfel M, Pérez de Prado A, et al. Three- and 6-month optical coherence tomographic surveillance following percutaneous coronary intervention with the Angiolite® drug-eluting stent: The ANCHOR study. *Catheter Cardiovasc Interv*. 2018;91:435-43.
13. Stone GW, Midei M, Newman W, Sanz M, Hermiller JB, Williams J, et al; SPIRIT III Investigators. Comparison of an everolimus-eluting stent and a paclitaxel-eluting stent in patients with coronary artery disease: a randomized trial. *JAMA*. 2008;299:1903-13.
14. Windecker S, Haude M, Neumann FJ, Stangl K, Witzensbichler B, Slagboom T, et al. Comparison of a novel biodegradable polymer sirolimus-eluting stent with a durable polymer everolimus-eluting stent: results of the randomized BIOFLOW-II trial. *Circ Cardiovasc Interv*. 2015;8:e001441.
15. Teeuwen K, van der Schaaf RJ, Adriaenssens T, Koolen JJ, Smits PC, Henriques JP, et al. Randomized multicenter trial investigating angiographic outcomes of hybrid sirolimus eluting stents with biodegradable polymer compared with everolimus-eluting stents with durable polymer in chronic total occlusions: The PRISON IV Trial. *JACC Cardiovasc Interv*. 2017;10:133-43.

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