



NeverTouch Gold-Tip Fiber
ENDOVENOUS LASER TREATMENT



The NeverTouch Gold-Tip Fiber for Less Pain and Bruising

The NeverTouch gold-tip laser fiber establishes a new standard for patient comfort, visibility and ease of use. The gold-tip fiber eliminates laser tip contact with the vein wall, which in turn minimizes perforations of the vein wall that typically result in pain and bruising as compared to traditional bare-tip fibers. The gold-tip also maximizes ultrasonic visibility, making it easier for physicians to use.




Echogenic



The Proprietary Science Behind NeverTouch Fibers

The innovation behind the NeverTouch fiber is a glass weld at the distal tip of a 600 μm fiber. The weld results in an effective fiber diameter of 905 μm and lowers the actual power density by 56 percent from that of a standard bare tip 600 μm fiber. The net effect is a homogeneous ablation with less focal charring of the vein wall than that which is seen with bare tip fibers. AngioDynamics' patent-pending technology is different from competitors with standard covered-tip fibers that have the same power density as a 600 μm fiber.



Procedure Kit Features and Benefits

Sheath-Lok system
Enables precise placement at the saphenofemoral junction by locking the fiber to the sheath, assuring the fiber stays a safe distance outside the sheath

Micro-access system in convenient "kangaroo pouch"
Select from a variety of micro-access kits, all packaged on the outside in a separate sterile pouch

4F TRE-Sheath hydrophilic-coated introducer (translucent, radiopaque, echogenic)
Echogenic for easy positioning and hydrophilic coating to reduce procedure time

Stainless steel braid
Improves pushability through tortuous anatomy and increases echogenicity

NeverTouch-FRS fiber
Fiber recognition system is compatible with all AngioDynamics' lasers

Sliding sheath gauge
Aids in marking treatment length and holding sheath in place during tumescent anesthesia

Available with a 21 or 19 gauge needle
New 19 gauge needle now available for physicians who prefer to insert the .035" wire directly through the needle

90 cm TRE-Sheath Introducer
Reach new lengths with the 90 cm NeverTouch procedure kit—ideal for gaining access below the knee

*Micro-Access System

AngioDynamics offers several micro-access systems to satisfy a variety of clinical needs. The chart below specifies the kit options and components that accompany the NeverTouch gold-tip fiber.

NeverTouch Fiber .018" System

The new .018" system eliminates the need for a 5F micro-access sheath and eliminates the need to exchange wires and size up, resulting in fewer procedure steps, and a faster procedure time.

	21 G Needle	19 G Needle	.018" Micro-access guidewire	Micro-introducer sheath with dilator	.035" long-access double-ended guidewire	.018" long-access guidewire	TRE-Sheath with dilator
Standard System	X		X	X	X		X
.018" System	X					X	X
19 G Micro-access System		X			X		X

AngioDynamics' VenaCure EVLT endovenous laser vein treatment offers patients a minimally-invasive choice for treating the source of their varicose veins and provides them with quicker recovery and a return to normal daily routines, as compared to surgical stripping. The VenaCure EVLT System includes a 1470 nm laser, 600 um or 400 um laser fiber procedure kits including accessories, marketing materials, support and more. For more information on these products, including our proprietary NeverTouch fiber technology, please visit www.VenaCure-EVLT.com.

NeverTouch Gold-Tip Fiber Procedure Kits

Description	UPN	Kit Length	Long-Access Guidewire	Needle
NeverTouch-FRS Kit	H787114030010	45 cm	.035"	21 gauge
NeverTouch-FRS Kit	H787114030020	65 cm	.035"	21 gauge
NeverTouch-FRS Kit	H787114040030	25 cm	.018"	21 gauge
NeverTouch-FRS Kit	H787114030040	90 cm	.035"	21 gauge
NeverTouch-FRS Kit	H787114030050	45 cm	.035"	19 gauge
NeverTouch-FRS Kit	H787114030060	65 cm	.035"	19 gauge
NeverTouch-FRS Kit	H787114030070	90 cm	.035"	19 gauge
NeverTouch-FRS Kit	H787114030120	45 cm	.018"	21 gauge
NeverTouch-FRS Kit	H787114030130	65 cm	.018"	21 gauge

TRE-Sheath Introducer Set

Description	UPN	Kit Length
TRE-Sheath Introducer Set	H787114033015	45 cm
TRE-Sheath Introducer Set	H787114033025	65 cm

Consult your local AngioDynamics representative for country specific availability.

IMPORTANT RISK INFORMATION

INDICATION FOR USE: The AngioDynamics, Inc. VenaCure EVLT NeverTouch Procedure Kits are indicated for endovascular coagulation of the Great Saphenous Vein (GSV) in patients with superficial vein reflux, for the treatment of varicose veins and varicosities associated with superficial reflux of the Great Saphenous Vein (GSV), and for the treatment of incompetence and reflux of superficial veins of the lower extremity. This product should be used only with lasers cleared for use in the treatment of varicose veins, varicosities with superficial reflux of the GSV, and in the treatment of incompetent refluxing veins in the superficial venous system in the lower limbs.

CONTRAINDICATIONS: Patients with thrombus in the vein segment to be treated. Patients with an aneurysmal section in the vein segment to be treated. Patients with peripheral artery disease as determined by an Ankle-Brachial Index <0.9. Patients with an inability to ambulate. Patients with deep vein thrombosis (DVT). Patients who are pregnant or breast feeding. Patients in general poor health. Other contraindications may be raised by the individual physician at the time of treatment.

PRECAUTION: Prior to and during use, avoid damaging the fiber by striking, stressing, or excessive bending. Do not coil the fiber tighter than a diameter of 20cm. Clinical safety and

effectiveness data is not available for other fiber tip designs and diameters. Prior to and during use, avoid bending the introducer sheath and dilator as this can cause kinks and damage. All persons in the treatment room must be wearing protective glasses with the proper rating for the wavelength being used.

WARNINGS: Treatment of a vein located close to the skin surface may result in skin burn. Paresthesia may occur from thermal damage to adjacent sensory nerves. Tissue not targeted for treatment must be protected from injury by direct and reflected laser energy. Use appropriate eye and protective wear for both patient and operating personnel. Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your sales representative. Inspect prior to use to verify that no damage has occurred during shipping. Storage instructions: Keep away from sunlight; Keep dry; Upper limit of temperature: 27 degrees Celsius. Reuse of single-use devices creates a potential risk of patient or user infections. Contamination of the device may lead to injury, illness, or death of the patient. Reprocessing may compromise the integrity of the device and/or lead to device failure, which, in turn, may result in patient injury, illness or death. VenaCure EVLT NeverTouch Procedure Kit is to be treated as contaminated biomedical waste subsequent to use. Additionally, the needle provided within the

kit is a sharp device. The used or unused devices, as well as the sharps, should be disposed of in accordance with hospital, administrative and/or local government policy for such items. Uncontaminated device packaging should be recycled if applicable, or disposed of as common waste in accordance with hospital, administrative and/or local government policy for such items. The fiber, guidewire and needle contain cobalt. Cobalt is classified as CMR 1B and is present in a concentration above 0.1% weight by weight.

POTENTIAL COMPLICATIONS: The potential complications include but are not limited to the following: vessel perforation, thrombosis, pulmonary embolism, phlebitis, hematoma, infection, skin pigmentation alteration, neovascularization, paresthesia due to thermal damage of adjacent sensory nerves, anesthetic tumescence, non-target

irradiation, hemorrhage, necrosis, air embolism, scarring, skin burns and pain.

Indications, contraindications, warnings and instructions for use can be found in the instructions for use supplied with each device. Observe all instructions prior to use. Failure to do so may result in patient complications.

CAUTION: Federal (USA) law restricts the sale of this device by or on the order of a physician.



USA > 14 Plaza Drive, Latham, NY 12110 > tel: 800-772-6446 or 518-798-1215 > fax: 518-798-1360
International > Haaksbergweg 75 (Margriettoren), 1101 BR, Amsterdam Z-O > The Netherlands
tel: +31 (0)20 753 2949 > fax: +31 (0)20 753 2939

www.angiodynamics.com | www.VenaCure-EVLT.com

AngioDynamics, the AngioDynamics logo, VenaCure EVLT, the VenaCure EVLT logo, NeverTouch, TRE-Sheath and Sheath-Lok are trademarks and/or registered trademarks of AngioDynamics, Inc., an affiliate or subsidiary.
© 2024 AngioDynamics, Inc. INV/BR/2140 Rev 01 02/2024

Distribuido por:

cardiva[®]

www.cardiva.com