

iVascular[®]

therapies for living

oceanus 35

Catéter balón de dilatación para ATP

Alcanzar & dilatar

Óptimo
desinflado

A detailed diagram of the Oceanus 35 catheter. The catheter is shown in a perspective view, highlighting its long, thin shaft and the balloon at the tip. The balloon is shown in a deflated state. Three callout lines with circular endpoints point to specific features: one points to the balloon, another to the double-lumen shaft, and a third to the conical tip. The background is a dark blue gradient.

Diseño de doble lumen
optimizado para una
excelente empujabilidad
y **navegabilidad**

Diseño de punta cónica
atraumática para un
óptimo cruce de
lesiones severas



www.ivascular.global

oceanus 35

características

- > Catéter coaxial (OTW)
- > Longitud útil del catéter: **80 cm, 140 cm y 200 cm**
- > Compatible con guía: **0.035"**
- > Perfil de punta: **0.036"**
- > Perfil de cruce: **0.057" a 0.083"**

Producto con Marca CE, certificado por Organismo Notificado 0318

Datos OCEANUS 35

Díámetro del balón (mm)	3mm	4mm	5mm	6mm	7mm	8mm	9mm	10mm	12mm
Presión nominal (atm)	7	7	7	6	7	7	6	6	6
Presión máx. recomendada (RBP) (atm)	16*	16	16	16	16	15	14	13	13
Presión de estallido promedio (ABP) (atm)	24	24	23	23	23	21	20	20	19

*20 atm: 3mm - L ≤ 60mm

Longitud útil del catéter	Balón Ø (mm)	Longitud del balón (mm)							
		20	40	60	80	120	150	200	
80 cm	3	BP PC35 080 300 020	BP PC35 080 300 040	BP PC35 080 300 060	BP PC35 080 300 080	BP PC35 080 300 120	BP PC35 080 300 150	BP PC35 080 300 200	
	4	BP PC35 080 400 020	BP PC35 080 400 040	BP PC35 080 400 060	BP PC35 080 400 080	BP PC35 080 400 120	BP PC35 080 400 150	BP PC35 080 400 200	
	5	BP PC35 080 500 020	BP PC35 080 500 040	BP PC35 080 500 060	BP PC35 080 500 080	BP PC35 080 500 120	BP PC35 080 500 150	BP PC35 080 500 200	
	6	BP PC35 080 600 020	BP PC35 080 600 040	BP PC35 080 600 060	BP PC35 080 600 080	BP PC35 080 600 120	BP PC35 080 600 150	BP PC35 080 600 200	
	7	BP PC35 080 700 020	BP PC35 080 700 040	BP PC35 080 700 060	BP PC35 080 700 080	BP PC35 080 700 120	-	-	
	8	BP PC35 080 800 020	BP PC35 080 800 040	BP PC35 080 800 060	BP PC35 080 800 080	-	-	-	
	9	BP PC35 080 900 020	BP PC35 080 900 040	BP PC35 080 900 060	BP PC35 080 900 080	-	-	-	
	10	BP PC35 080 100 020	BP PC35 080 100 040	BP PC35 080 100 060	-	-	-	-	
	12	BP PC35 080 120 020	BP PC35 080 120 040	BP PC35 080 120 060	-	-	-	-	
140 cm	3	BP PC35 140 300 020	BP PC35 140 300 040	BP PC35 140 300 060	BP PC35 140 300 080	BP PC35 140 300 120	BP PC35 140 300 150	BP PC35 140 300 200	
	4	BP PC35 140 400 020	BP PC35 140 400 040	BP PC35 140 400 060	BP PC35 140 400 080	BP PC35 140 400 120	BP PC35 140 400 150	BP PC35 140 400 200	
	5	BP PC35 140 500 020	BP PC35 140 500 040	BP PC35 140 500 060	BP PC35 140 500 080	BP PC35 140 500 120	BP PC35 140 500 150	BP PC35 140 500 200	
	6	BP PC35 140 600 020	BP PC35 140 600 040	BP PC35 140 600 060	BP PC35 140 600 080	BP PC35 140 600 120	BP PC35 140 600 150	BP PC35 140 600 200	
	7	BP PC35 140 700 020	BP PC35 140 700 040	BP PC35 140 700 060	BP PC35 140 700 080	BP PC35 140 700 120	-	-	
	8	BP PC35 140 800 020	BP PC35 140 800 040	BP PC35 140 800 060	BP PC35 140 800 080	-	-	-	
	9	BP PC35 140 900 020	BP PC35 140 900 040	BP PC35 140 900 060	BP PC35 140 900 080	-	-	-	
	10	BP PC35 140 100 020	BP PC35 140 100 040	BP PC35 140 100 060	-	-	-	-	
	12	BP PC35 140 120 020	BP PC35 140 120 040	BP PC35 140 120 060	-	-	-	-	
NEW 200 cm	3	-	-	-	-	-	-	-	
	4	-	-	-	-	-	-	-	
	5	-	-	-	-	-	-	-	
	6	BP PC35 200 600 020	-	-	-	-	-	-	
	7	BP PC35 200 700 020	-	-	-	-	-	-	
	8	BP PC35 200 800 020	-	-	-	-	-	-	
	9	BP PC35 200 900 020	-	-	-	-	-	-	
	10	BP PC35 200 100 020	-	-	-	-	-	-	
	12	BP PC35 200 120 020	-	-	-	-	-	-	

● Compatibilidad con introductor: **5F** ● Compatibilidad con introductor: **6F**

Distribuido por:



Fabricado por:

Life Vascular Devices Biotech, S.L.
www.ivascular.global
info@ivascular.global